sigvaris

Compreboot® Standard Foot Compreboot® Plus Foot

Gebrauchsanweisung	2	Petunjuk Penggunaan	30
Instructions for Use	3	使用指南	32
Istruzioni per l'uso	5	사용 설명서	33
Mode d'emploi	6	使用方法	34
Brugsanvisning	8	คำแนะนํ้าในการใช้งาน	35
Käyttöohje	9	Пайдалану жөнiндегі нұсқаулық	36
Gebruiksaanwijzing	11	उपयोग के लिए निर्देश	38
Bruksanvisning	12	Arahan Penggunaan	39
Instrucciones de uso	13	40	تعليمات الاستخدام
Instruções de utilização	15	42	הנחיות לשימוש
Инструкция по применению	16		
Instrukcja użytkowania	18		
Návod na použitie	20		
Navodila za uporabo	21		
Οδηγίες χρήσης	22		
Kullanım talimatı	24		
Инструкции за употреба	25		
Instrucțiuni de utilizare	27		
Naudojimo instrukcija	29		

Produktbeschreibung

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English

Product description

Compreboot® Standard Foot and Compreboot® Plus Foot provide up to 40 mmHg of foot compression and are made from a soft, supple material.

Intended purpose

These inelastic wraps are compression therapy devices used on the foot. Compression therapy improves venous and lymphatic return, which reduces venous and lymphatic edema. It is used in the initial decongestion phase of lymphatic edema and pronounced venous edema as an alternative to wrapping with bandages. This medical device is intended for multiple use and is not sterile.

Effect of medical compression

Compression products apply a carefully controlled pressure on tissue and veins. This improves the blood return in the veins and increases the lymphatic flow. At the same time, it improves microcirculation and therefore the supply to the tissue.

Target group

The intended target group includes patients who have been prescribed a compression product for individual compression therapy by the responsible physician or trained health professional. The target group also includes user groups from health care professions who apply the products to patients.

Indications and contraindications

Indications

- Edema caused by chronic venous disease (CVD) due to obstruction or valvular insufficiency; treatment and/or prevention
- CVD with venous edema, skin changes, open and closed ulcers
- Lymphedema
- Lipedema
- Edema due to immobilization
- Non-specific edema or edema after surgery

Contraindications

- Advanced peripheral arterial occlusive disease PAOD (ABPI <0.6, ankle blood pressure <60 mmHg, toe blood pressure <30 mmHg, or TcPO₂ <20 mmHg instep)
- Severe cardiac insufficiency (NYHA III-IV)

- Septic phlebitis
- Phlegmasia coerulea dolens
- Severe sensory disturbances of the extremity
- Advanced peripheral neuropathy (e.g. in diabetes mellitus)
- Compression of an existing arterial bypass (epifascial)

Relative contraindications that require a particular balancing of benefits and risks:

- Allergy against compression material
- Inflammatory processes (compression products to reduce inflammation, pain or edema in leg erysipelas or cellulitis should only be used in combination with antibacterial treatment)
- PAOD with ABPI between 0.6 and 0.9 and/or ankle pressure between 60 and 90 mmHg
- Skin diseases with massive exudation
- Primary chronic polyarthritis

Risks and side effects

Inelastic wraps can cause skin necrosis and pressure damage to the peripheral nerves if handled improperly. Itching, scaling and signs of inflammation may occur under the compression aids in cases of sensitive skin, so adequate skin care is advised. Regular medical check-ups are required to maintain successful therapy. If acute pain or skin irritation occurs while wearing, immediately consult your physician. In this case, the compression aid should no longer be worn until otherwise prescribed by a physician.

Please note: Consult your physician. The precautions and contraindications listed here are guidelines only and are no substitute for medical advice.

Special instructions

- Protect the compression product from ointments, oils and greases to preserve its service life.
- Compression products must not be worn over open wounds. Wounds should be covered with appropriate dressings.
- Failure to comply with this instruction may reduce the safety and effectiveness of the product.
- All serious incidents that occur when using the product should be reported to the manufacturer and the competent supervisory authority.
- Compression is increased when liners with compression pressure are used under a wrap.
- It should be ensured that the pressure of the compression product decreases (degressive pressure gradient) from distal to proximal (away from the body to close to the body).

Application

Applying the Compreboot® Standard Foot/ Compreboot® Plus Foot

1. Make sure that the hook fasteners are open.
2. Place your heel in the heel area of the wrap.
3. Close the hook fasteners around your foot.
4. The hook fasteners may need to be readjusted several times to achieve a comfortable fit. Make sure that stockings,

socks etc. which are worn under the wrap do not have any creases

5. Readjust the product as needed to ensure proper compression.

Removing the Compreboot® Standard Foot/ Compreboot® Plus Foot

To remove the product, loosen the hook fasteners starting at the ankle, roll their ends backwards and fasten them there. Once all the hook fasteners have been secured, remove the product from the foot.

Storage, wearing period and service life

Unless otherwise prescribed by your physician, the compression product can be worn during the day and overnight. However, remove the product once a day to clean the skin and care for any wounds. The product should then be reapplied according to the instructions for use. How often the product needs to be removed for skin checks, hygiene and/or wound care may vary according to the individual needs of patients as well as the prescriber's recommendations. A liner or stocking should be worn under every compression product. Wearing and washing a compression product every day may cause the medical pressure gradient and elasticity of your wrap to diminish. The shelf life date is indicated on the packaging using the symbol of an hourglass. We guarantee maximum effect for the subsequent six months of wearing.

Care and storage

Follow the instructions for care on the packaging and/or the textile label. Close all hook fasteners before washing. Wash the product in a mesh laundry bag on a gentle cycle using a mild detergent. Do not use bleach or fabric softener. Dry the product on low/no heat in a dryer or air dry (air drying is recommended to extend wrap life). Store the product at room temperature in a dry place away from light. Keep away from direct heat.

Disposal

There are no special requirements for disposing of the medical compression product. It can be disposed of as residual or household waste. Please comply with local regulations for the disposal of textiles.

List of materials

Compreboot® Standard Foot

Neoprene
Polyamide
Elastane
Made without natural rubber latex

Compreboot® Plus Foot

Polyamide
Polyurethane foam
Neoprene
Polyester
Elastane
Made without natural rubber latex

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SIGVARIS and SIGVARIS GROUP are registered trademarks of SIGVARIS AG,
© 2023 Copyright by SIGVARIS AG, CH-9014 St. Gallen/Switzerland
SIGVARIS S.r.l., Via Galileo Galilei 2 / 2, 31057 Silea (TV), **Italy**
SIGVARIS S.A., ul. Powidzka 50, Gutowo Małe, 62-300, Września, **Poland**
SIGVARIS (Australia) Pty. Ltd., 8/53 Stanley Street, West Melbourne VIC 3003, **Australia**
SIGVARIS (SHANGHAI) TRADING CO LTD, Pudong New District, Unit 207, Building 1,