sigvaris

Compreflex® Standard Arm

Gebrauchsanweisung	2	Petunjuk Penggunaan	30
Instructions for Use	3	使用指南	31
Istruzioni per l'uso	5	사용 설명서	32
Mode d'emploi	6	使用方法	33
Brugsanvisning	8	คำแนะนํ้าในการใช้งาน	35
Käyttöohje	9	Пайдалану жөнiндегi	
Gebruiksaanwijzing	10	нұсқаулық	36
Bruksanvisning	12	उपयोग के लिए निर्देश	37
Instrucciones de uso	13	Arahan Penggunaan	38
Instruções de utilização	15	40	تعليمات الاستخدام
Инструкция по применению	16	41	הנחיות לשימוש
Instrukcja użytkowania	18		
Návod na použitie	20		
Navodila za uporabo	21		
Οδηγίες χρήσης	22		
Kullanım talimatı	24		
Инструкции за употреба	25		
Instrucțiuni de utilizare	27		
Naudojimo instrukcija	28		

Produktbeschreibung

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English

Product description

Compreflex® Standard Arm is an inelastic product that provides a mild to medium compression (≤ 40 mmHg). The product is combined with a pair of Cotton Liners to protect the skin.

Intended purpose

These inelastic wraps are compression therapy devices used on the arm. Compression therapy improves venous and lymphatic return, which reduces venous and lymphatic edema. It is used in the initial decongestion phase of lymphatic edema and pronounced venous edema as an alternative to wrapping with bandages. This medical device is intended for multiple use and is not sterile.

Effect of medical compression

Compression products apply a carefully controlled pressure on tissue and veins. This improves the blood return in the veins and increases the lymphatic flow. At the same time, it improves microcirculation and therefore the supply to the tissue.

Target group

The intended target group includes patients who have been prescribed a compression product for individual compression therapy by the responsible physician or trained health professional. The target group also includes user groups from health care professions who apply the products to patients.

Indications and contraindications

Indications

- Edema
- Lymphedema
- Lipedema

Contraindications

- Advanced peripheral arterial occlusive disease PAD (radial arterial pressure < 60 mmHg, finger arterial pressure < 30 mmHg, or TcPO₂ dorsum of the hand < 20 mmHg)
- Severe cardiac insufficiency (NYHA III-IV)
- Septic phlebitis
- Severe sensory disturbances of the extremity

- Advanced peripheral neuropathy (e.g. in diabetes mellitus)
- Compression of an existing arterial bypass (epifascial)

Relative contraindications that require a particular balancing of benefits and risks:

- Allergy against compression material
- Inflammatory processes (compression products to reduce inflammation, pain or edema due to erysipelas or cellulitis should only be used in combination with antibacterial treatment)
- Mild PAD (reduced radial arterial pressure in lateral comparison, but >60 mmHg radial arterial pressure)
- Skin diseases with massive exudation
- Blue or white coloration of the fingers, numbness of the hands
- Primary chronic polyarthritis

Risks and side effects

Inelastic wraps can cause skin necrosis and pressure damage to the peripheral nerves if handled improperly. Itching, scaling and signs of inflammation may occur under the compression aids in cases of sensitive skin, so adequate skin care is advised. Regular medical check-ups are required to maintain successful therapy. If acute pain or skin irritation occurs while wearing, immediately consult your physician. In this case, the compression aid should no longer be worn until otherwise prescribed by a physician.

Please note: Consult your physician. The precautions and contraindications listed here are guidelines only and are no substitute for medical advice.

Special instructions

- Protect the compression product from ointments, oils and greases to preserve its service life.
- Compression products must not be worn over open wounds. Wounds should be covered with appropriate dressings.
- Failure to comply with this instruction may reduce the safety and effectiveness of the product.
- All serious incidents that occur when using the product should be reported to the manufacturer and the competent supervisory authority.
- It should be ensured that the pressure of the compression product decreases (degressive pressure gradient) from distal to proximal (away from the body to close to the body).
- Make sure that the liner does not crease.

Application

Applying the Compreflex® Standard Arm

1. First, please remove all jewelry from the hand, wrist and arm. Pull the provided liner on over your arm and up to the shoulder. Make sure that the knitted section is evenly distributed over the entire length
2. Before applying, unfold each hook fastener on the product to its full length, roll the ends back and secure them there so that they cannot get tangled. Pull the product carefully over the arm, but no further than

the liner, and align the elbow seam to the elbow.

3. Close all hook fasteners starting from the wrist without applying compression. The hook fasteners should overlap as little as possible without leaving any gap between them.
4. Then tighten all hook fasteners to a comfortable compression level.
5. Finally, attach butterfly strap to elbow crease area.

Removing the Compreflex® Standard Arm

First remove the butterfly strap. Then undo every hook fastener starting from the wrist. Roll their ends backwards and secure them there so that they cannot get tangled. Remove the product from the arm and take off the liner.

Storage, wearing period and service life

Unless otherwise prescribed by your physician, the compression product can be worn during the day and overnight. However, remove the product once a day to clean the skin and care for any wounds. The product should then be reapplied according to the instructions for use. How often the product needs to be removed for skin checks, hygiene and/or wound care may vary according to the individual needs of patients and the prescriber's recommendations. A liner should be worn under every compression product. Wearing and washing a compression product every day may cause the medical pressure gradient and elasticity of your wrap to diminish. The shelf life date is indicated on the packaging using the symbol of an hourglass. We guarantee maximum effect for the subsequent six months of wearing.

Care and storage

Follow the instructions for care on the packaging and/or the textile label. Close all hook fasteners before washing. Wash the product in a mesh laundry bag on a gentle cycle using a mild detergent. Do not use bleach or fabric softener. Dry the product on low/no heat in a dryer or air dry (air drying is recommended to extend wrap life). Store the product at room temperature in a dry place away from light. Keep away from direct heat. Wash liners before first use.

Disposal

There are no special requirements for disposing of the medical compression product. It can be disposed of as residual or household waste. Please comply with local regulations for the disposal of textiles.

List of materials

Polyamide
 Polyurethane foam
 Elastane
 Cotton
 Neoprene
 Made without natural rubber latex

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SIGVARIS and SIGVARIS GROUP are registered trademarks of SIGVARIS AG,
© 2023 Copyright by SIGVARIS AG, CH-9014 St. Gallen/Switzerland
SIGVARIS S.r.l., Via Galileo Galilei 2 / 2, 31057 Silea (TV), **Italy**
SIGVARIS S.A., ul. Powidzka 50, Gutowo Małe, 62-300, Września, **Poland**
SIGVARIS (Australia) Pty. Ltd., 8/53 Stanley Street, West Melbourne VIC 3003, **Australia**
SIGVARIS (SHANGHAI) TRADING CO LTD, Pudong New District, Unit 207, Building 1,