sigvaris

Compreflex® Standard Calf and Foot

Gebrauchsanweisung	2	Petunjuk Penggunaan	36
Instructions for Use	4	使用指南	38
Istruzioni per l'uso	5	사용 설명서	39
Mode d'emploi	7	使用方法	40
Brugsanvisning	9	คำแนะนํ้าในการใช้ขา	41
Käyttöohje	11	Пайдалану жөнiндегi	
Gebruiksaanwijzing	12	нұсқаулық	43
Bruksanvisning	14	उपयोग के लिए निर्देश	45
Instrucciones de uso	16	Arahan Penggunaan	46
Instruções de utilização	18	48	تعليمات الاستخدام
Инструкция по применению	20	49	הנחיות לשימוש
Instrukcja użytkowania	21		
Návod na použitie	23		
Navodila za uporabo	25		
Οδηγίες χρήσης	27		
Kullanım talimatı	29		
Инструкции за употреба	30		
Instrucțiuni de utilizare	32		
Naudojimo instrukcija	34		

Produktbeschreibung

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Product description

The Compreflex® Standard Calf and Foot is an inelastic two-component product that provides precise compression ranges of 21–46 mmHg in the calf region thanks to the use of the Accutabs. The element provides up to 40 mmHg of compression for the foot. The product is combined with a pair of Basic Liners to protect the skin.

Intended purpose

These inelastic wraps are compression therapy devices used on the calf and foot. Compression therapy improves venous and lymphatic return, which reduces venous and lymphatic edema. It is used in the initial decongestion phase of lymphatic edema and pronounced venous edema as an alternative to wrapping with bandages. This medical device is intended for multiple use and is not sterile.

Effect of medical compression

Compression products apply a carefully controlled pressure on tissue and veins. This improves the blood return in the veins and increases the lymphatic flow. At the same time, it improves microcirculation and therefore the supply to the tissue.

Target group

The intended target group includes patients who have been prescribed a compression product for individual compression therapy by the responsible physician or trained health professional. The target group also includes user groups from health care professions who apply the products to patients.

Indications and contraindications

Indications

- Edema caused by chronic venous disease (CVD) due to obstruction or valvular insufficiency; treatment and/or prevention
- CVD with venous edema, skin changes, open and closed ulcers
- Lymphedema
- Lipedema
- Edema due to immobilization
- Non-specific edema or edema after surgery

Contraindications

- Advanced peripheral arterial occlusive disease PAOD (ABPI <0.6, ankle blood pressure <60 mmHg, toe blood pressure <30 mmHg, or TcPO₂ <20 mmHg instep)
- Severe cardiac insufficiency (NYHA III–IV)
- Septic phlebitis
- Phlegmasia coerulea dolens
- Severe sensory disturbances of the extremity
- Advanced peripheral neuropathy (e.g. in diabetes mellitus)
- Compression of an existing arterial bypass (epifascial)

Relative contraindications that require a particular balancing of benefits and risks:

- Allergy against compression material
- Inflammatory processes (compression products to reduce inflammation, pain or edema

in leg erysipelas or cellulitis should only be used in combination with antibacterial treatment)

- PAOD with ABPI between 0.6 and 0.9 and/or ankle pressure between 60 and 90 mmHg
- Skin diseases with massive exudation
- Primary chronic polyarthritis

Risks and side effects

Inelastic wraps can cause skin necrosis and pressure damage to the peripheral nerves if handled improperly. Itching, scaling and signs of inflammation may occur under the compression aids in cases of sensitive skin, so adequate skin care is advised. Regular medical check-ups are required to maintain successful therapy. If acute pain or skin irritation occurs while wearing, immediately consult your physician. In this case, the compression aid should no longer be worn until otherwise prescribed by a physician. **Please note:** Consult your physician. The precautions and contraindications listed here are guidelines only and are no substitute for medical advice.

Special instructions

- Protect the compression product from ointments, oils and greases to preserve its service life.
- Compression products must not be worn over open wounds. Wounds should be covered with appropriate dressings.
- Failure to comply with this instruction may reduce the safety and effectiveness of the product.
- All serious incidents that occur when using the product should be reported to the manufacturer and the competent supervisory authority.
- Compression is increased when liners with compression pressure are used under a wrap.
- It should be ensured that the pressure of the compression product decreases (degressive pressure gradient) from distal to proximal (away from the body to close to the body).
- Make sure that the liner does not crease.

Application

Applying the Compreflex® Standard Calf and Foot

Pull the liner (included) evenly over your heel and distribute the knitted section evenly over the leg. The optimum fit of the liner is achieved when it finishes 2 finger widths below the back of the knee.

Foot element

1. Make sure that the hook fasteners are open.
2. Place your heel in the heel area of the wrap.
3. Close the hook fasteners around your foot.
4. The hook fasteners may need to be readjusted several times to achieve a comfortable fit.
5. Readjust the product as needed to ensure proper compression.

Calf element

1. Before applying, unfold each hook fastener on the product to its full length, roll the

ends back and secure them there so that they cannot get tangled. Slide the foot into the opening between the elastic insert and the product and pull the product up to just below the knee, but no further than the liner. The textile label is located on the top at the back of the leg and faces outward.

- Carefully position the two lowest hook fasteners over the leg and fasten them without applying compression. The hook fasteners should overlap as little as possible without leaving any gap between them. Press the ends of the fasteners firmly to the sides. Continue positioning the remaining hook fasteners in the same way and fix them in place.
- Apply the Accutabs: your prescribing physician should tell you which compression range is best for you. Apply the Accutabs accordingly and then close the hook fasteners as directed (refer to the Accutab® instructions for use).

Removing the Compreflex® Standard Calf and Foot

To remove the product, loosen the hook fasteners starting at the ankle, roll their ends backwards and fasten them there. Once all the hook fasteners have been secured, remove the product from the foot. Also loosen each of the calf compression aid's hook fasteners, roll their ends backwards and secure them there so that they cannot get tangled. Remove the product from the leg and take off the liner.

Storage, wearing period and service life

Unless otherwise prescribed by your physician, the compression product can be worn during the day and overnight. However, remove the product once a day to clean the skin and care for any wounds. The product should then be reapplied according to the instructions for use. How often the product needs to be removed for skin checks, hygiene and/or wound care may vary according to the individual needs of patients as well as the prescriber's recommendations. A liner should be worn under every compression product. Wearing and washing a compression product every day may cause the medical pressure gradient and elasticity of your wrap to diminish. The shelf life date is indicated on the packaging using the symbol of an hourglass. We guarantee maximum effect for the subsequent six months of wearing.

Care and storage

Follow the instructions for care on the packaging and/or the textile label. Close all hook fasteners before washing. Wash the product in a mesh laundry bag on a gentle cycle using a mild detergent. Do not use bleach or fabric softener. Dry the product on low/no heat in a dryer or air dry (air drying is recommended to extend wrap life). Store the product at room temperature in a dry place away from light. Keep away from direct heat. Wash liners before first use.

Disposal

There are no special requirements for disposing of the medical compression product. It can be disposed of as residual or household waste. Please comply with local regulations for the disposal of textiles.

List of materials

Neoprene
Polyamide
Elastane
Polyurethane foam
Cotton
Polyester
Made without natural rubber latex

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SIGVARIS and SIGVARIS GROUP are registered trademarks of SIGVARIS AG,
© 2023 Copyright by SIGVARIS AG, CH-9014 St. Gallen/Switzerland
SIGVARIS S.r.l., Via Galileo Galilei 2 / 2, 31057 Silea (TV), **Italy**
SIGVARIS S.A., ul. Powidzka 50, Gutowo Małe, 62-300, Września, **Poland**
SIGVARIS (Australia) Pty. Ltd., 8/53 Stanley Street, West Melbourne VIC 3003, **Australia**
SIGVARIS (SHANGHAI) TRADING CO LTD, Pudong New District, Unit 207, Building 1,